



Pantex

Division of Bio-Analysis, Inc.

Providing Diagnostic Technology for a Better Tomorrow.

CERTIFICATE OF ANALYSIS:

ESTRIOL (E3) Rabbit Polyclonal Antibody

Catalog Number: E3-373
Lot Number: 373
Clonality: Polyclonal
Host: Rabbit
Immunogen: Estriol-6-carboxymethyloxime -BSA
Testing Method: EIA
Supplied as: Neat serum antibody

| Compound | |
|---|----------------------------------|
| <i>C-18 Steroids:</i> | <i>% Cross-reactivity</i> |
| Estriol | 100.00 |
| Estradiol-17 β | 0.0550 |
| Estriol-3-glucuronide | 2.2917 |
| Estriol-3-SO ₄ | 0.8462 |
| Estriol-16-glucuronide | 0.0275 |
| Estradiol-17 β -3-SO ₄ | 0.0149 |
| 17-Epiestriol | 0.0611 |
| 16,17-Epiestriol | 0.1196 |
| 16-Hydroxyestrone | 0.9167 |
| 16-Ketoestradiol | 0.2391 |
| 15-Hydroxyestriol | 32.3529 |
| Estriol-17-SO ₄ | 0.1122 |

Rev: 2