

Certificate of Analysis

Product:	Anti-Testosterone Polyclonal Antibody
Catalog Number:	AB190
Lot No.:	336F
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	Testosterone-19-CME-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
C-19 Steroids:	% Cross-reactivity
Testosterone	100.00
5α-Dihydrotestosterone	1.780
Androsterone	0.001
Dehydroepiandrosterone	0.002
Dehydroepiandrosterone-SO ₄	< 0.001
Androstenedione	1.140
Androstenediol	0.880
C-21 Steroids:	% Cross-reactivity
Progesterone	0.020
17α-OH-Progesterone	0.050
Pregnenolone	< 0.001
17-OH-Pregnenolone	< 0.001
Desoxycorticosterone	0.020
Corticosterone	0.002
Cortisol	0.001
11-Desoxycortisol	0.006
Aldosterone	
C-18 Steroids:	% Cross-reactivity
Estradiol-17β	<0.001
Estradiol-17α	0.005
Estrone	< 0.001
Estriol	< 0.001

Rev: 3