

Certificate of Analysis

Product:	Anti-Corticosterone Polyclonal Antibody
Catalog Number:	AB120
Lot No.:	277P
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	Corticosterone-3-hemisuccinate-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
<i>C-21 Steroids:</i>	<i>% Cross-reactivity</i>
Corticosterone	100.00
Aldosterone	< 0.001
Cortisol	0.056
11-Desoxycortisol	0.061
Desoxycorticosterone	1.970
Pregnenolone	< 0.001
Progesterone	0.001
17-OH-Pregnenolone	< 0.001
17-OH-Progesterone	< 0.001
<i>C-19 Steroids:</i>	<i>% Cross-reactivity</i>
Androstanediol	< 0.001
Androstenedione	< 0.001
Androsterone	< 0.001
Dehydroepiandrosterone	< 0.001
Dehydroepiandrosterone-SO ₄	0.010
5α-Dihydrotestosterone	< 0.001
Testosterone	< 0.001
<i>C-18 Steroids:</i>	<i>% Cross-reactivity</i>
Estradiol-17β	< 0.001
Estradiol-17α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 4