

Certificate of Analysis

Product:	Anti-DHEA-S Polyclonal Antibody
Catalog Number:	AB140
Lot No.:	201P
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	DHEA-3-S04-7-CMO-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
C-19 Steroids:	% Cross-reactivity
Dehydroepiandrosterone-SO ₄	100.00
Dehydroepiandrosterone	< 0.0003
Testosterone	0.0007
5α-Dihydrotestosterone	0.0029
Androstenedione	0.0483
C-21 Steroids:	% Cross-reactivity
Progesterone	< 0.0003
17α-OH-Progesterone	0.0023
Pregnenolone	< 0.0003
17-OH-Pregnenolone	0.0014
17-OH-Pregnenolone-SO ₄	< 0.0003
Desoxycorticosterone	0.0003
Corticosterone	< 0.0003
Cortisol	< 0.0003
11-Desoxycortisol	0.0133
C-18 Steroids:	% Cross-reactivity
Estradiol-17β	< 0.0003
Estradiol-17α	< 0.0003
Estrone	< 0.0003
Estriol	< 0.0003

Rev: 3

Certificate of Analysis

Product:	Anti-DHEA-S Polyclonal Antibody
Catalog Number:	AB140
Lot No.:	204P
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	DHEA-3-S04-7-CMO-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
<i>C-19 Steroids:</i>	<i>% Cross-reactivity</i>
Dehydroepiandrosterone-SO ₄	100.00
Dehydroepiandrosterone	< 0.001
Testosterone	< 0.001
5α-Dihydrotestosterone	< 0.001
Androstenedione	< 0.001
Androsterone	< 0.001
<i>C-21 Steroids:</i>	<i>% Cross-reactivity</i>
Progesterone	< 0.001
17α-OH-Progesterone	< 0.001
Pregnenolone	< 0.001
17-OH-Pregnenolone	< 0.001
17-OH-Pregnenolone-SO ₄	< 0.001
Desoxycorticosterone	< 0.001
Corticosterone	< 0.001
Cortisol	< 0.001
11-Desoxycortisol	< 0.001
Aldosterone	< 0.001
<i>C-18 Steroids:</i>	<i>% Cross-reactivity</i>
Estradiol-17β	< 0.001
Estradiol-17α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 3