

Certificate of Analysis

Product:	Anti-Triiodothyronine (T3) Monoclonal Antibody
Catalog Number:	MAB290
Lot/No.:	7
Clonality:	Monoclonal
Host Species:	Rabbit
Immunogen:	T3-HSA
Source:	Recombinant DNA (sequenced from original Hybridoma)
Testing Method:	EIA
Supplied as:	mg of IgG Antibody (purified)
Buffer:	PBS with 40% Glycerol
Storage:	≤ -15°C
Concentration:	1mg/1mL
Expiration Date:	As specified on label
Specificity:	For the quantitative determination of Total T3, Free T3 (FT3) and Neonatal T3

Compound:	Concentration ng/mL	Cross-reactivity
L-Triiodothyronine	-	1.0
L-Thyroxine	1,000	0.0126
3,5 Diiodo-DL-Tyrosine	10,000	≤0.0001
3-Iodo-L-Tyrosine	10,000	≤0.0001

Rev: 1