



Pantex

Division of Bio-Analysis, Inc.

Providing Diagnostic Technology for a Better Tomorrow.

CERTIFICATE OF ANALYSIS:

CORTICOSTERONE Rabbit Polyclonal Antibody

Catalog Number: CORTICO-277P
Lot Number: 277P
Clonality: Polyclonal
Host: Rabbit
Immunogen: Corticosterone-3-hemisuccinate-BSA
Testing Method: EIA
Supplied as: Neat serum antibody

Compound	
C-21 Steroids:	% Cross-reactivity
Corticosterone	100.00
Aldosterone	< 0.001
Cortisol	0.056
11-Desoxycortisol	0.061
Desoxycorticosterone	1.970
Pregnenolone	< 0.001
Progesterone	0.001
17-OH-Pregnenolone	< 0.001
17-OH-Progesterone	< 0.001
C-19 Steroids:	% Cross-reactivity
Androstenediol	< 0.001
Androstenedione	< 0.001
Androsterone	< 0.001
Dehydroepiandrosterone	< 0.001
Dehydroepiandrosterone-SO ₄	0.010
5 α -Dihydrotestosterone	< 0.001
Testosterone	< 0.001
C-18 Steroids:	% Cross-reactivity
Estradiol-17 β	< 0.001
Estradiol-17 α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 3