

Certificate of Analysis

Product: Anti-Estriol (E3) Polyclonal Antibody
Catalog Number: AB171
Lot No.: 373
Clonality: Polyclonal
Host Species: Rabbit
Immunogen: Estriol-6-carboxymethyloxime -BSA
Testing Method: EIA
Supplied as: Neat serum antibody
Storage: ≤ -15°C
Expiration Date: As specified on label

Compound	
<i>C-18 Steroids:</i>	<i>% Cross-reactivity</i>
Estriol	100.00
Estradiol-17β	0.0550
Estriol-3-glucuronide	2.2917
Estriol-3-SO ₄	0.8462
Estriol-16-glucuronide	0.0275
Estradiol-17β-3-SO ₄	0.0149
17-Epiestriol	0.0611
16,17-Epiestriol	0.1196
16-Hydroxyestrone	0.9167
16-Ketoestradiol	0.2391
15-Hydroxyestriol	32.3529
Estriol-17-SO ₄	0.1122

Rev: 3