

Certificate of Analysis

Product:	Anti-Aldosterone Monoclonal Antibody
Catalog Number:	MAB330
Lot/No.:	134-1
Clonality:	Monoclonal
Host Species:	Rabbit
Immunogen:	Aldosterone-3-CMO-BSA
Source:	Recombinant DNA (sequenced from original Hybridoma)
Testing Method:	EIA
Supplied as:	mg of IgG Antibody (purified)
Buffer:	PBS with 40% Glycerol
Storage:	≤ -15°C
Concentration:	1mg/1mL
Expiration Date:	As specified on label

Compound:	% Cross-reactivity
Aldosterone	100.00
Cortisol	< 0.001
Cortisone	0.012
Corticosterone	0.010
Progesterone	< 0.001
DHEA-SO ₄	0.016
Estradiol-17β	0.008
Estriol	0.008

Rev: 3