



Pantex

Division of Bio-Analysis, Inc.

Providing Diagnostic Technology for a Better Tomorrow.

CERTIFICATE OF ANALYSIS:

ALDOSTERONE Rabbit Monoclonal Antibody

Catalog Number: ALDOSTERONE-134-1
Lot Number: 134-1
Clonality: **Monoclonal**
Immunogen: Aldosterone-3-CMO-BSA
Source: **Recombinant DNA** (sequenced from original Hybridoma)
Testing Method: EIA
Supplied as: µg or mg of IgG Antibody (purified)
Buffer: PBS with 40% Glycerol
Storage: ≤ -15°C
Concentration: As specified on label
Expiration Date: As specified on label

Compound:	% Cross-reactivity
Aldosterone	100.00
Cortisol	< 0.001
Cortisone	0.012
Corticosterone	0.010
Progesterone	< 0.001
DHEA-SO ₄	0.016
Estradiol-17β	0.008
Estriol	0.008

Rev: 2