

Certificate of Analysis

Product:	Anti-Cortisol Polyclonal Antibody
Catalog Number:	AB110
Lot No.:	24P
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	Cortisol-3-0-carboxymethyloxime-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
<i>C-21 Steroids:</i>	<i>% Cross-reactivity</i>
Cortisol	100.00
17α-OH Progesterone	0.645
Pregnenolone	< 0.001
17-OH-Pregnenolone	< 0.001
Progesterone	0.0222
Desoxycorticosterone	0.0544
11-Desoxycortisol	3.425
Dexamethasone	0.0103
Cortisone	2.083
Corticosterone	0.660
Aldosterone	< 0.001
<i>C-19 Steroids:</i>	<i>% Cross-reactivity</i>
Androstenedione	< 0.001
Testosterone	< 0.001
5α-Dihydrotestosterone	< 0.001
Dehydroepiandrosterone SO4	< 0.001
Androstanedione	< 0.001
<i>C-18 Steroids:</i>	<i>% Cross-reactivity</i>
Estradiol-17β	< 0.001
Estradiol-17α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 3

Certificate of Analysis

Product:	Anti-Cortisol Polyclonal Antibody
Catalog Number:	AB110
Lot No.:	425A
Clonality:	Polyclonal
Host Species:	Rabbit
Immunogen:	Cortisol-3-0-carboxymethyloxime-BSA
Testing Method:	EIA
Supplied as:	Neat serum antibody
Storage:	≤ -15°C
Expiration Date:	As specified on label

Compound	
<i>C-21 Steroids:</i>	<i>% Cross-reactivity</i>
Cortisol	100.00
17α-OH Progesterone	0.8966
Pregnenolone	< 0.001
17-OH-Pregnenolone	< 0.001
Progesterone	0.0361
Desoxycorticosterone	0.2826
11-Desoxycortisol	7.222
Dexamethasone	0.0057
Cortisone	2.6000
Corticosterone	1.9697
Aldosterone	< 0.001
<i>C-19 Steroids:</i>	<i>% Cross-reactivity</i>
Androstenedione	< 0.001
Testosterone	< 0.001
5α-Dihydrotestosterone	< 0.001
Dehydroepiandrosterone-SO ₄	< 0.001
Androstenedione	< 0.001
<i>C-18 Steroids:</i>	<i>% Cross-reactivity</i>
Estradiol-17β	< 0.001
Estradiol-17α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 3