



Pantex

Division of Bio-Analysis, Inc.

Providing Diagnostic Technology for a Better Tomorrow.

CERTIFICATE OF ANALYSIS:
ESTRONE (E1) Rabbit Monoclonal Antibody

Catalog Number:	ESTRONE (E1) – 94-2
Lot Number:	94-2
Clonality:	Monoclonal
Immunogen:	Estrone-6-CMO-BSA
Source:	Recombinant DNA (sequenced from original Hybridoma)
Testing Method:	EIA
Supplied as:	µg or mg of IgG Antibody (purified)
Buffer:	PBS with 40% Glycerol
Storage:	≤ -15°C
Concentration:	As specified on label
Expiration Date:	As specified on label

Compounds	
C-18 Steroids:	% Cross-reactivity
Estrone	100.00
Estradiol-17β	0.8282
Estradiol-17α	0.2107
Estriol	0.0347
C-19 Steroids:	% Cross-reactivity
Dehydroepiandrosterone-SO4	0.0017
C-21 Steroids:	% Cross-reactivity
Progesterone	0.0055

Rev: 2