



Pantex

Division of Bio-Analysis, Inc.

Providing Diagnostic Technology for a Better Tomorrow.

CERTIFICATE OF ANALYSIS:

TESTOSTERONE Rabbit Polyclonal Antibody

Catalog Number: TESTO-336F
Lot Number: 336F
Clonality: Polyclonal
Host: Rabbit
Immunogen: Testosterone-19-CME-BSA
Testing Method: EIA
Supplied as: Neat serum antibody

Compound	
C-19 Steroids:	% Cross-reactivity
Testosterone	100.00
5 α -Dihydrotestosterone	1.780
Androsterone	0.001
Dehydroepiandrosterone	0.002
Dehydroepiandrosterone-SO ₄	< 0.001
Androstenedione	1.140
Androstenediol	0.880
C-21 Steroids:	% Cross-reactivity
Progesterone	0.020
17 α -OH-Progesterone	0.050
Pregnenolone	< 0.001
17-OH-Pregnenolone	< 0.001
Desoxycorticosterone	0.020
Corticosterone	0.002
Cortisol	0.001
11-Desoxycortisol	0.006
Aldosterone	
C-18 Steroids:	% Cross-reactivity
Estradiol-17 β	<0.001
Estradiol-17 α	0.005
Estrone	< 0.001
Estriol	< 0.001

Rev: 2