

Certificate of Analysis

Product: Anti-DHEA Polyclonal Antibody
Catalog Number: AB130
Lot No.: 147P
Clonality: Polyclonal
Host Species: Rabbit
Immunogen: DHEA-15-Carboxyethylthioether-BSA
Testing Method: EIA
Supplied as: Neat serum antibody
Storage: ≤ -15°C
Expiration Date: As specified on label

Compound	
C-19 Steroids:	% Cross-reactivity
Dehydroepiandrosterone	100.00
Dehydroepiandrosterone-SO ₄	0.0037
Testosterone	0.002
5α-Dihydrotestosterone	0.007
Androstenedione	0.056
C-21 Steroids:	% Cross-reactivity
Progesterone	0.230
17α-OH-Progesterone	0.0004
Pregnenolone	0.013
17-OH-Pregnenolone	0.072
Desoxycorticosterone	0.052
Corticosterone	0.004
Cortisol	0.0007
11-Desoxycortisol	0.012
Aldosterone	0.0003
C-18 Steroids:	% Cross-reactivity
Estradiol-17β	< 0.001
Estradiol-17α	< 0.001
Estrone	< 0.001
Estriol	< 0.001

Rev: 3